

All-Ceramic Cementation

BRUX – Full Contour Zirconia

- Ceramir Crown and Bridge Cement
- Hybrid glass ionomer cement (RelyX Luting Cement)
- 3M ESPE; GCFuji Plus, GC America)
- Resin cements for short or over-tapered preparations (RelyX Unicem, 3M ESPE; Panavia F2.0, Kuraray)

e.Max – Lithium Disilicate

- Below is a list of recommended cements:
Vivaglass (Ivoclar Vivadent)
GC Fuji (GC America)
Ketec Cem (3M ESPE)
Panavia F (J. Morita)
C&B Metabond
Variolink (Ivoclar Vivadent)

Note: Vitremer Is Not Recommended For Any All-ceramic Restoration

Layered Zirconia

- Panavia 21
- Glass Ionomer Cement (Fuji GC)
- Resin Ionomer Cement (RelyX by 3M ESPE)